

Quick guide for Physicians on SleepLink

S.M.A. is now called SleepLink. While SleepLink has a different look and feel, the basic functions and underlying data are the same. You can use your S.M.A. login and password to enter SleepLink. This quick guide will help you get started with using SleepLink. If you have any questions please call 1-800-334-5085.

The screenshot shows the 'Physician - Studies' page in a web browser. The URL is <https://sma.sleepmed.md/PhysStudyList.aspx>. The page has a navigation menu with 'My Profile', 'Quality Review', 'HST Manager', 'Physician', 'HST Physician', and 'Management Reports'. A 'Logout' button is in the top right. The main heading is 'Physician - Studies' with the SleepMed logo. Search filters include 'Lab' (Pasadena), 'From' (4/12/2013), 'To' (6/12/2013), 'Physician' (a dropdown menu), and 'Status' (2 items checked). There are 'Search' and 'Export to Excel' buttons. A table below shows search results with columns for Patient Name, Lab, IV #, Status, Study Type, Night Tech, Scorer, Study Date, Scored Date, and Interp. Date. Two callouts provide additional information: one points to the Physician dropdown stating 'Physicians can now be picked from a dropdown box only showing physicians you are associated with.' and another points to the Status dropdown stating 'A Search can include multiple statuses, physicians, or even labs. For example, both Scored and In Interpretation statuses are selected here.' At the bottom, there are 'Retrieve Selected Studies' and 'Retrieve Selected Reports' buttons, a disclaimer about disabled checkboxes, and a copyright notice for SleepMed Administrator.

Physicians can now be picked from a dropdown box only showing physicians you are associated with.

A Search can include multiple statuses, physicians, or even labs. For example, both Scored and In Interpretation statuses are selected here.

If you wish to view studies that have already been signed, please select the interpreted status and click search.

	Patient Name	Lab	IV #	Status	Study Type	Night Tech	Scorer	Study Date	Scored Date	Interp. Date
<input type="checkbox"/>	[REDACTED]	Pasadena	01898465	Scored	Split	murakami, keith	Tran, Hoa	06/03/2013	06/04/2013	
<input type="checkbox"/>	[REDACTED]	Pasadena	01913503	In Interpretation	Split	Kerwat, Mitchell	Tran, Hoa	06/07/2013	06/10/2013	
<input type="checkbox"/>	[REDACTED]	Pasadena	01910789	Scored	Split	Meave, Lupe	Priddy, Michael	06/10/2013	06/11/2013	
<input type="checkbox"/>	[REDACTED]	Pasadena	01913683	Scored	Split		Priddy, Michael	06/07/2013	06/11/2013	
<input type="checkbox"/>	[REDACTED]	Pasadena	01896115	Scored	Split	Bui, Duy	Priddy, Michael	06/10/2013	06/11/2013	

Retrieve Selected Studies Retrieve Selected Reports

* Disabled checkboxes indicate that the study data cannot be downloaded. Only studies that have been scored can be downloaded.

SleepMed Administrator Copyright © 2005 SleepMed incorporated

Physician - Studies



Filters: Location: Pasadena | From: 4/12/2013 | To: 6/12/2013 | Physician: | Status: 3 items checked | Scorer: Type in a Scorer | Night Tech: Type in a tech

Search | Export to Excel

If you wish to view studies that have already been signed, please select the interpreted status and click search.

Page size: 20 | 292 items in 15 pages

	Patient Name	Lab	IV #	Status	Study Type	Night Tech	Scorer	Scored Date	Interp. Date
<input type="checkbox"/>	[Redacted]	Pasadena	01841530	Interpreted	Split	Meave, Lupe	[Redacted]	04/17/2013	04/17/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01885438	Interpreted		Tran, Hoa	[Redacted]	05/20/2013	06/08/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01890272	Interpreted	PSG	Cazares, Jesus	[Redacted]	06/03/2013	06/04/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01855378	Interpreted		Meave, Lupe	Pickett, Debra	05/03/2013	05/06/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01886857	Interpreted		Dang, Tai	Tran, Hoa	05/09/2013	05/10/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01868542	Interpreted		murakami, keith	Pickett, Debra	04/22/2013	04/23/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01865142	Interpreted		Bui, Duy	Pickett, Debra	04/16/2013	04/17/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01878575	Interpreted	Split	Meave, Lupe	Reseigh, Steven	05/07/2013	05/08/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01892362	Interpreted	Split		Priddy, Michael	05/21/2013	05/23/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01868024	Interpreted	Split	Meave, Lupe	Tran, Hoa	06/03/2013	06/05/2013
<input type="checkbox"/>	[Redacted]	Pasadena	01891550	Interpreted	Split	Kerwat, Mitchell	Priddy, Michael	05/22/2013	05/23/2013

Click on any column header to sort all patients in ABC or 123 order.

Control the number of studies that appear on each page here.

After a Search, results can be imported into Microsoft Excel.

The screenshot shows a web browser window with the URL https://sma.sleepmed.md/sma_temppro/PhysStudyFormReport.aspx?StudyID=08231507-4434-4785-af7e-466924d75a7b&ID=&ControlName=Physician&IV=01898465&AcctNo=519754&St. The browser's address bar shows the page title "Physician - Studies".

The browser's menu bar includes "File", "Edit", "View", "Favorites", "Tools", and "Help". The left sidebar contains a "Reports" section with a list of documents, including "SPLIT Dr. Saraf 6-3-2013 71506 PSW.doc", and an "Images" section with a list of items such as "I.V.", "Sleep-CMN 5/20/13", "1st F2F - 05.20.2013", "TECH NOTES - 5/28/13", "TECH NOTES - 06/3/13", "Sleep-BILLING, DOS- 6/3/13", and "Request for CMN-DME 06.15.2013".

The main content area displays a report with a "Menu" dropdown menu open, listing options: "Save report to server", "Retrieve study", "Send study", "Sign report", "Print report", "Delete report", "Send study back for re-scoring", "Save report locally", and "Attach report".

The report text includes a "PROCEDURE" section: "An all-night comprehensive sleep study was performed in which the following medical parameters were recorded using a Healthdyne model Alice 3 and Alice 5 computerized polygraph, left central (C3), right central (C4), left occipital (O1) and right occipital (O2) electroencephalogram, left and right electrooculogram, electrocardiogram, submental & anterior tibialis electromyogram, nasal/oral airflow, oxygen saturation (pulse oximetry), chest & abdominal effort belts, sonogram (snoring) and body position monitor. The study raw data was manually reviewed and interpreted by Pratap Saraf, M.D., Diplomate of the American Board of Internal Medicine. The recording started at 6/3/2013 at 9:06:00 PM, 6/4/2013 at 5:06:00 AM."

The report also includes a "DEFINITION" section with a bulleted list of criteria for respiratory events.

Callouts provide instructions: "Click on Menu to Sign Report and all other familiar options." (pointing to the Menu dropdown), "Access Tech notes and Insurance information. Click on one to view it." (pointing to the Reports sidebar), and "Click on the arrow here to collapse the left pane and maximize the report or image currently being viewed." (pointing to a red arrow icon on the left side of the report content area).

Click on Menu to Sign Report and all other familiar options.

Access Tech notes and Insurance information. Click on one to view it.

Click on the arrow here to collapse the left pane and maximize the report or image currently being viewed.

